

HWP-A Organic Wheat Hydrolysate

Product Code: 6304



QMS Certification Services

Organic for use in high value nutritional applications

HWP-A is an enzymatically hydrolyzed organic wheat gluten powder. This material is filtered for improved flowability and filterability in a variety of applications. This material is completely soluble when used at concentrations as high as 35%. This version of HWP-A is Organic, Animal-free, Kosher and HALAL.

SPECIFICATIONS:

Typical Analysis:

Moisture	<6.0%
Total Nitrogen (dry basis)	≥ 11.2%
Amino Nitrogen (dry basis)	≥ 2.8%
pH (10% solution)	6.0-7.3
Standard Plate Count	<5,000 cfu/g
Coliform	<10 cfu/g
Yeast and Mold	<50 cfu/g
Salmonella	Negative/25g

Packaging:


25 Kg Polyethylene bag in a plastic drum

Storage Requirements:

This product should be kept in a cool, dry and ventilated place. The product should be used as soon as possible after the bag has been opened. Any partially used bags must be sealed properly to prevent moisture absorption by product.

Shelf Life:

The shelf life of the product is 3 years when stored at <30° C.



952.936.3600

5400 Opportunity Court,
Suite 120
Minnetonka, MN 55343

www.Nu-TekBioSciences.com